

Product datasheet for TP300640

DUSP4 (NM_001394) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dual specificity phosphatase 4 (DUSP4), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200640 protein sequence Red =Cloning site Green =Tags(s) MVTMEELREMDCSVLKRLMNRDENGSGAGSGSHGTLGLPSGGKCLLLDCRPFLAHSAGYILGSVNVRC N TIVRRRAKGSVSLEQILPAEEVRRRLRSLGLYSAVIVYDERSPRAESLREDSTVSLVVQALRRNAERTDI CLLKGGYERFSSEYPEFCSTKALAAIPPPVPPSATEPLDLGCSSCGTPLHDQGGPVEILPFLYLGSAYH AARRDMLDALGITALLNVSSDCPNHFEGHYQYKCIPVEDNHNKADISSWFMEIAIEYIDAVKDCRGRVLVHC QAGISRSATICLAYLMMKKRVRLEEAFFVKQRRSIISPNFSFMGQLLQFESQVLATSCAAEAASPSGPL RERGKTPATPTSQVFVSPVSVGVHSAPSSLPYLHSPITTSPSC TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	42.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: [NP_001385](#)

Locus ID: 1846

UniProt ID: [Q13115](#)

RefSeq Size: 5625

Cytogenetics: 8p12

RefSeq ORF: 1182

Synonyms: HVH2; MKP-2; MKP2; TYP

Summary: The protein encoded by this gene is a member of the dual specificity protein phosphatase subfamily. These phosphatases inactivate their target kinases by dephosphorylating both the phosphoserine/threonine and phosphotyrosine residues. They negatively regulate members of the mitogen-activated protein (MAP) kinase superfamily (MAPK/ERK, SAPK/JNK, p38), which are associated with cellular proliferation and differentiation. Different members of the family of dual specificity phosphatases show distinct substrate specificities for various MAP kinases, different tissue distribution and subcellular localization, and different modes of inducibility of their expression by extracellular stimuli. This gene product inactivates ERK1, ERK2 and JNK, is expressed in a variety of tissues, and is localized in the nucleus. Two alternatively spliced transcript variants, encoding distinct isoforms, have been observed for this gene. In addition, multiple polyadenylation sites have been reported. [provided by RefSeq, Jul 2008]

Protein Families: Phosphatase

Protein Pathways: MAPK signaling pathway

Product images:



Coomassie blue staining of purified DUSP4 protein (Cat# TP300640). The protein was produced from HEK293T cells transfected with DUSP4 cDNA clone (Cat# [RC200640]) using MegaTran 2.0 (Cat# [TT210002]).