

Product datasheet for **TP300648L**

Aldehyde dehydrogenase 10 (ALDH3A2) (NM_001031806) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aldehyde dehydrogenase 3 family, member A2 (ALDH3A2), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200648 protein sequence Red =Cloning site Green =Tags(s)
	<p>MELEVRRVRQAFLSGRSRPLRFRLQQLALRRMVQEREKDLTAIAADLCKSEFNVYSQEVITVLGEIDF MLENLPEWVTAKPVKKNVLTMLDEAYIQPQPLGVVLIIGAWNYPFVLTIQPLIGAIAAGNAVIKPSSELS ENTAKILAKLLPQYLDQDLYIVINGGVEETELLKQRFDHIFYTGNTAVGKIVMEAAAKHLTPVTLELGG KSPCYIDKDCDLDIVCRITWGYMNCGQTCIAPDYILCEASLQNQIVWKIKETVKEFYGENIKESPDYE RIINLRHFKRILSLLEGQKIAFGGETDEATRYIAPTDLTDVDPKTKVMQEEIFGPILPIVPVKNVDEAIN FINEREKPLALYVFSHNHKLKRMIDETSSGGVTGNDVIMHFTLNSFPFGGSSGMGAYHGKHSFDTFE HQRPCLLKSLKREGANKLRYPPNSQSKVDWGFLLKRFNKEKLGLLLLTFLGIVAAVLVKKYQAVLRRK ALLIFLVHRLRWSSKQR</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	57.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001026976
Locus ID:	224
UniProt ID:	P51648
RefSeq Size:	3823
Cytogenetics:	17p11.2
RefSeq ORF:	1524
Synonyms:	ALDH10; FALDH; SLS
Summary:	Aldehyde dehydrogenase isozymes are thought to play a major role in the detoxification of aldehydes generated by alcohol metabolism and lipid peroxidation. This gene product catalyzes the oxidation of long-chain aliphatic aldehydes to fatty acid. Mutations in the gene cause Sjogren-Larsson syndrome. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Arginine and proline metabolism, Ascorbate and aldarate metabolism, beta-Alanine metabolism, Butanoate metabolism, Fatty acid metabolism, Glycerolipid metabolism, Glycolysis / Gluconeogenesis, Histidine metabolism, Limonene and pinene degradation, Lysine degradation, Metabolic pathways, Propanoate metabolism, Pyruvate metabolism, Tryptophan metabolism, Valine, leucine and isoleucine degradation

Product images:



Coomassie blue staining of purified ALDH3A2 protein (Cat# [TP300648]). The protein was produced from HEK293T cells transfected with ALDH3A2 cDNA clone (Cat# [RC200648]) using MegaTran 2.0 (Cat# [TT210002]).