

Product datasheet for **TP300658**

Syntaxin 3 (STX3) (NM_004177) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human syntaxin 3 (STX3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200658 protein sequence Red =Cloning site Green =Tags(s)

MKDRLEQLKAKQLTQDDDDAVEIAIDNTAFMDEFFSEIEETRLNIDKISEHVVEEAKKLYSIILSAPIPE
PKTKDDLEQLTTEIKKRANNVRNKLKSMEKHIEEDEVRSADLRIRKSQHSVLSRKFVEVMTKYNEAQVD
FRERSKGRIQRQLEITGKKTDEELEEMLESGNPAIFTSGIIDSQISKQALSEIEGRHKDIVERLESSIKE
LHDMFMDIAMLVENQGEMLDNIELNVMHTVDHVEKARDETKKAVKYQSQARKKLIIIVLVVLLGILAL
IIGLSVGLN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	33 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004168
Locus ID:	6809



[View online »](#)

UniProt ID: [Q13277](#), [Q53YE2](#)

RefSeq Size: 6498

Cytogenetics: 11q12.1

RefSeq ORF: 867


Synonyms: STX3A

Summary: The gene is a member of the syntaxin family. The encoded protein is targeted to the apical membrane of epithelial cells where it forms clusters and is important in establishing and maintaining polarity necessary for protein trafficking involving vesicle fusion and exocytosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010]

Protein Families: Druggable Genome, Stem cell - Pluripotency, Transmembrane

Protein Pathways: SNARE interactions in vesicular transport

Product images:



Coomassie blue staining of purified STX3 protein (Cat# TP300658). The protein was produced from HEK293T cells transfected with STX3 cDNA clone (Cat# [RC200658]) using MegaTran 2.0 (Cat# [TT210002]).