

Product datasheet for **TP300673M**

PFKP (NM_002627) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphofructokinase, platelet (PFKP), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200673 protein sequence Red =Cloning site Green =Tags(s)

MDADDSRAPKGLRKFLEHLSGAGKAIGVLTSGGDAQGMNAAVRAWVRMGIYVVGAKVYFIYEGYQGMVDG
GSNIAEADWESVSSILQVGGTIIGSARCQAFRTREGRLKAACNLLQRGITNLCVIGGDGSLTGANLFRKE
WSSLLEELARNGQIDKEAVQKYAYLNVVGMVGSIDNDFCGTDMTIGTDSALHRIIEVVDAIMTTAQSHQR
TFVLEVMGRHCGYLALVSALACGADWWFLPESPPEEGWEEQMCVKLSENRRARKRLNIIIVAEGAIDTQN
KPITSEKIKELVWTQLGYDTRVTILGHVQRGGTPSAFDRILASRMGVEAVIALLEATPDTACVSLNGN
HAVRLPLMECVQMTQDVQKAMDERRFQDAVRLRGRSFAGNLNTYKRLAIKLPDDQIPKTNCNVAVINVGA
PAAGMNAAVRSVAVRVGIADGHRMLAIYDGFDFGFAKGGQIKIIGWTDVGGWTGQGGGSLGTRKRVLP GKYLEE
IATQMRTHSINALLIIGGFEAYLGLLELSAAREKHEEFCVPMVMVPATVSNVPGSDFSIGADTALNTIT
DTCDRIKQSASGTRRRVFIETMGGYCGYLANMGGLAAGADAAYIFEFPDIRDLQSNVEHLTEKMKTTI
QRGLVLRNESCENYTTDFIYQLYSEEGKGVDFCRKNVLGHMQGGGAPSPFDRNFGTKISARAMEWITAK
LKEARGRGKFTTDDSIKRVIFQPVAELKKQTD FEHRIPKEQWWLKLRLPMKILAKYKASYDV
SDSGQLEHVQPWSV

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	85.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Cell treatment (PMID: 28467929)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.




[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002618
Locus ID:	5214
UniProt ID:	Q01813
RefSeq Size:	2657
Cytogenetics:	10p15.2
RefSeq ORF:	2352
Synonyms:	ATP-PFK; PFK-C; PFK-P; PFKF
Summary:	This gene encodes a member of the phosphofructokinase A protein family. The encoded enzyme is the platelet-specific isoform of phosphofructokinase and plays a key role in glycolysis regulation. This gene may play a role in metabolic reprogramming in some cancers, including clear cell renal cell carcinomas, and cancer of the bladder, breast, and lung. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]
Protein Families:	Druggable Genome
Protein Pathways:	Fructose and mannose metabolism, Galactose metabolism, Glycolysis / Gluconeogenesis, Metabolic pathways, Pentose phosphate pathway

Product images:



PFKP promotes the proliferation of human pulmonary artery smooth muscle cells (HPASMCs). HPASMCs were serum-starved for 48 hours and stimulated with recombinant PFKP protein (OriGene [TP300673]) (100 ng/mL) for 72 hours. Fixed cells stained with crystal violet were counted to assess cell proliferation. Figure cited from Cell Metab, PMID: 28467929



Coomassie blue staining of purified PFKP protein (Cat# [TP300673]). The protein was produced from HEK293T cells transfected with PFKP cDNA clone (Cat# [RC200673]) using MegaTran 2.0 (Cat# [TT210002]).