

## Product datasheet for **TP300709L**

### MEA1 (NM\_014623) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human male-enhanced antigen 1 (MEA1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200709 protein sequence <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MGPERHLSGAPARMATVVLGGDTMGPERIFPNQTEELGHQGPSEGTGDWSSEEPEEEQEETGSGPAGYS Y QPLNQDPEQEEVELAPVGDGDVADIQDRIQALGLHLPDPPLSEDEDEEGATALNNHSSIPMDPEHVEL VKRTMAGVSLPAPGVPAWAREISDAQWEDVWQKALQARQASPAWK  <span style="color: red;">TR</span> <span style="color: green;">TRPLEQKLISEEDLAANDILDYKDDDDKV</span>
Tag:	C-Myc/DDK
Predicted MW:	19.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_055438</a></u>
Locus ID:	4201
UniProt ID:	<u><a href="#">Q16626</a></u>


[View online »](#)

RefSeq Size:	868
Cytogenetics:	6p21.1
RefSeq ORF:	555
Synonyms:	HYS; MEA
Summary:	May play an important role in spermatogenesis and/or testis development.[UniProtKB/Swiss-Prot Function]

### Product images:



Coomassie blue staining of purified MEA1 protein (Cat# [TP300709]). The protein was produced from HEK293T cells transfected with MEA1 cDNA clone (Cat# [RC200709]) using MegaTran 2.0 (Cat# [TT210002]).