

## Product datasheet for **TP300713L**

### **FKBP52 (FKBP4) (NM\_002014) Human Recombinant Protein**

#### **Product data:**

**Product Type:** Recombinant Proteins  
**Description:** Recombinant protein of human FK506 binding protein 4, 59kDa (FKBP4), 1 mg  
**Species:** Human  
**Expression Host:** HEK293T  
**Expression cDNA Clone or AA Sequence:** >RC200713 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MTAEEMKATESGAQSAPLPMEGVDISPKQDEGVLKVIKREGTGTEMPMIGDRVHVHYTGWLLDGTKFDSS  
LDRKDKFSFDLKGGEVIKAWDIAIATMKVGEVCHITCKPEYAYGSAGSPPKIPP NATLVFEVELFEFKGE  
DLTEEDGGIIRRIQTRGEGYAKPNEGAIVEVALEGYKDKLFDQRELFEIGEGENLDLPYGLERAIQR  
MEKGESIVYLKPSYAFGSVGKEKFQIPPNAELKYELHLKSFKAKESWEMNSEEKLEQSTIVKERGTVY  
FKEGKYKQALLQYKKIVSWLEYESSFSNEEAQKAQALRLASHLNLAMCHLKLQAFSAAIESCNKALELDS  
NNEKGLFRRGEAHLAVNDFELARADFQKVLQLYPNNKAAKTQLAVCQQRIRRLAREKKLYANMFERLAE  
EENKAKAEASSGDHPTDTEMKEEQKSNTAGSQSQVETEA

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK  
**Predicted MW:** 51.6 kDa  
**Concentration:** >0.05 µg/µL as determined by microplate BCA method  
**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining  
**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol  
**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.  
**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.  
**Storage:** Store at -80°C.  
**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.  
**RefSeq:** [NP\\_002005](#)




[View online >](#)

Locus ID:	2288
UniProt ID:	<a href="#">Q02790</a>
RefSeq Size:	3757
Cytogenetics:	12p13.33
RefSeq ORF:	1377
Synonyms:	FKBP51; FKBP52; FKBP59; HBI; Hsp56; p52; PPIase

**Summary:** The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds to the immunosuppressants FK506 and rapamycin. It has high structural and functional similarity to FK506-binding protein 1A (FKBP1A), but unlike FKBP1A, this protein does not have immunosuppressant activity when complexed with FK506. It interacts with interferon regulatory factor-4 and plays an important role in immunoregulatory gene expression in B and T lymphocytes. This encoded protein is known to associate with phytanoyl-CoA alpha-hydroxylase. It can also associate with two heat shock proteins (hsp90 and hsp70) and thus may play a role in the intracellular trafficking of hetero-oligomeric forms of the steroid hormone receptors. This protein correlates strongly with adeno-associated virus type 2 vectors (AAV) resulting in a significant increase in AAV-mediated transgene expression in human cell lines. Thus this encoded protein is thought to have important implications for the optimal use of AAV vectors in human gene therapy. The human genome contains several non-transcribed pseudogenes similar to this gene. [provided by RefSeq, Sep 2008]

**Protein Families:** Druggable Genome

### Product images:



Coomassie blue staining of purified FKBP4 protein (Cat# [TP300713]). The protein was produced from HEK293T cells transfected with FKBP4 cDNA clone (Cat# [RC200713]) using MegaTran 2.0 (Cat# [TT210002]).