

Product datasheet for TP300726

RRM1 (NM_001033) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ribonucleotide reductase M1 (RRM1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>Peptide sequence encoded by RC200726 Blue=ORF Red=Cloning site Green=Tag(s)

MHVIKRDGRQERVMFDKITSRIQKLCYGLNMDFVDPAQITMKVIQGLYSGVTTVELDTLAAETAATLTT
KHPDYAILAARIAVSNLHKETKKVFSVDMEDLYNYINPHNGKHSPMVAKSTLDIVLANKDRLNSAIIYD
RDFSYNFYFGFKTLERSYLLKINGKVAERPQHMLMRVSVGIHKEDIDAAIETYNLLSERWFTHASPTLFN
AGTNRQLSSCFLLSMKDDSIIEGIYDTLKQCALISKSAGGIGVAVSCIRATGSYIAGTNGNSNGLVPML
RVYNNTARYVDQGGNKRPGAFAIYLEPWHLDFEFDLKKNKTGKEEQRRARDLFFALWIPDLFMKRVETN
QDWSLMCPNECPGLDEWGWEEFEKLYASYEKQGRVRKVVKAQQLWYAIIESQTETGTPYMLYKDCSNRK
SNQQNLGTIKCSNLCTEIVEYTSKDEVAVCNLASLALNMYVTSEHTYDFKLAEVTKVVRNLNKIIDI
NYPVPEACLSNKRHRPIGIGVQGLADAFILMRYPFESAEAQLLNKQIFETIYYGALEASCDLAKEQGP
YETIEGSPVSKGILQYDMWNVPTDLWDWKVLKEKIAKYGIRNSLLIAPMPTASTAQILGNNESIEPYT
SNIYTRRVLSGEFQIVNPHLLKDLTERGLWHEEMKNQIIACNGSIQSIPEDDLKQLYKTVWEISQKT
VLKMAAERGAFIDQSQSLNIIHIAEPNYGKLTSMHFYGWKQGLKTMGYLLRTPAANPIQFTLNKEKLDK
KEKVSKEEEEKERNTAAMVCSLENRDECLMCGS
TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Recombinant protein using RC200726 also available, [TP300726](#)

Tag:	C-Myc/DDK
Predicted MW:	89.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	Binding assay (PMID: 29765556)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online >](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001024](#)

Locus ID: 6240

UniProt ID: [P23921](#)

RefSeq Size: 3117

Cytogenetics: 11p15.4

RefSeq ORF: 2376


Synonyms: R1; RIR1; RR1

Summary: This gene encodes the large and catalytic subunit of ribonucleotide reductase, an enzyme essential for the conversion of ribonucleotides into deoxyribonucleotides. A pool of available deoxyribonucleotides is important for DNA replication during S phase of the cell cycle as well as multiple DNA repair processes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015]


Protein Families: Druggable Genome

Protein Pathways: Glutathione metabolism, Metabolic pathways, Purine metabolism, Pyrimidine metabolism

Product images:



Afatinib, a covalent inhibitor of ErbB family, modifies RRM1 and RRM2. Two ug recombinant RRM1 (OriGene TP300726) and/or 1 ug RRM2 in the presence of 12.5 uM afatinib were incubated for one hour at 37 C. The reaction product was examined by Coomassie blue G-250 staining and immunoblotting with anti-afatinib antiserum. Figure cited from Oncotarget, PMID: 29765556



Coomassie blue staining of purified RRM1 protein (Cat# TP300726). The protein was produced from HEK293T cells transfected with RRM1 cDNA clone (Cat# [RC200726]) using MegaTran 2.0 (Cat# [TT210002]).