

## Product datasheet for **TP300747L**

### **MAD3 (MXD3) (NM\_031300) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human MAX dimerization protein 3 (MXD3), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200747 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MEPLASNIQVLLQAAEFLEERREREAHGYASLCPHRSPGPIHRRKKRPPQAPGAQDSGRSVHNELEKRRR AQLKRCLERLKQQMPLGADCARYTTLSLLRRARMHIQKLEDQEQRARQLKERLRSKQQSLQRQLEQLRGL AGAAERERLRADSLDSSGLSSERSDSDQEELEVDVESLVFGGEALLRGFVAGQEHSYSHGGGAWL  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	23.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_112590</a></u>
Locus ID:	83463
UniProt ID:	<u><a href="#">Q9BW11</a></u>
RefSeq Size:	1483



[View online »](#)

Cytogenetics: 5q35.3


RefSeq ORF: 618

Synonyms: BHLHC13; MAD3; MYX

**Summary:** This gene encodes a member of the Myc superfamily of basic helix-loop-helix leucine zipper transcriptional regulators. The encoded protein forms a heterodimer with the cofactor MAX which binds specific E-box DNA motifs in the promoters of target genes and regulates their transcription. Disruption of the MAX-MXD3 complex is associated with uncontrolled cell proliferation and tumorigenesis. Transcript variants of this gene encoding different isoforms have been described.[provided by RefSeq, Dec 2008]

**Protein Families:** Druggable Genome, Transcription Factors

### Product images:



Coomassie blue staining of purified MXD3 protein (Cat# [TP300747]). The protein was produced from HEK293T cells transfected with MXD3 cDNA clone (Cat# [RC200747]) using MegaTran 2.0 (Cat# [TT210002]).