

Product datasheet for TP300757L

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SOX2 (NM 003106) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human SRY (sex determining region Y)-box 2 (SOX2), 1 mg

Species: Human Expression

Host:

HEK293T

Expression >RC200757 protein sequence cDNA Clone or Red=Cloning site Green=Tags(s) AA Sequence:

MYNMMETELKPPGPQQTSGGGGGNSTAAAAGGNQKNSPDRVKRPMNAFMVWSRGQRRKMAQENPKMHNSE

ISKRLGAEWKLLSETEKRPFIDEAKRLRALHMKEHPDYKYRPRRKTKTLMKKDKYTLPGGLLAPGGNSMA

SGVGVGAGLGAGVNQRMDSYAHMNGWSNGSYSMMQDQLGYPQHPGLNAHGAAQMQPMHRYDVSALQYNSM TSSQTYMNGSPTYSMSYSQQGTPGMALGSMGSVVKSEASSSPPVVTSSSHSRAPCQAGDLRDMISMYLPG

AEVPEPAAPSRLHMSQHYQSGPVPGTAINGTLPLSHM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 34.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Bioactivity: EMSA assay (PMID: 25892221)

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss

of protein during the filtration process.

Storage: Store at -80°C.

Stable for 12 months from the date of receipt of the product under proper storage and handling Stability:

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 003097





SOX2 (NM_003106) Human Recombinant Protein - TP300757L

Locus ID: 6657

UniProt ID: <u>P48431</u>, <u>A0A0U3FYV6</u>

RefSeq Size: 2520 Cytogenetics: 3q26.33

RefSeq ORF: 951

Synonyms: ANOP3; MCOPS3

Summary: This intronless gene encodes a member of the SRY-related HMG-box (SOX) family of transcription


factors involved in the regulation of embryonic development and in the determination of cell fate. The product of this gene is required for stem-cell maintenance in the central nervous system, and also regulates gene expression in the stomach. Mutations in this gene have been associated with optic nerve hypoplasia and with syndromic microphthalmia, a severe form of structural eye malformation. This gene lies within an intron of another gene called SOX2 overlapping transcript (SOX2OT). [provided by

RefSeq, Jul 2008]

Protein Adult stem cells, Cancer stem cells, Embryonic stem cells, ES Cell Differentiation/IPS, Induced

Families: pluripotent stem cells, Transcription Factors

Product images:



Coomassie blue staining of purified SOX2 protein (Cat# [TP300757]). The protein was produced from HEK293T cells transfected with SOX2 cDNA clone (Cat# [RC200757]) using MegaTran 2.0 (Cat# [TT210002]).