

## Product datasheet for TP300769M

### FAM110A (NM\_001042353) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 110, member A (FAM110A), transcript variant 3, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200769 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MPVHTLSPGAPSAPALPCRLRTRVPGYLLRGPADGGARKPSAVERLEADKAKYVKSLSLVANTRQEPVQPL LSKQPLFSPETRRTVLTPSRRALPGPCRRPQLDLIDLSSLDLDCDSPVSPAEASRTPGRAEGAGRPPPAT PPRPPPSTSAVRRVDVRLPASPARPCSPGPAAASSPARPPGLQRSKSDLSERFSSRAAADLERFFNFCG LDPEEARGLGVAHLARASSDIVSLAGPSAGPGSSEGGCSRSSVTVEERARERVPYGVSVVERNARVIKW LYGLRQARESPAEEG
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	31.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001035812</a>
Locus ID:	83541



[View online »](#)

UniProt ID: [Q9BQ89](#)  
RefSeq Size: 1827  
Cytogenetics: 20p13  
RefSeq ORF: 885  
Synonyms: bA371L19.3; C20orf55; F10

**Product images:**

Coomassie blue staining of purified FAM110A protein (Cat# [TP300769]). The protein was produced from HEK293T cells transfected with FAM110A cDNA clone (Cat# [RC200769]) using MegaTran 2.0 (Cat# [TT210002]).