

Product datasheet for **TP300780**

ZNF747 (NM_023931) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human zinc finger protein 747 (ZNF747), 20 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC200780 protein sequence
Red=Cloning site **Green**=Tags(s)

MTDPSLGLTVPMAPPLAPLPPRDPNGAGSEWRKPGAVSFADVAVYFSREEWGCLRPAQRALYRDVMRETY
GHLGALGESPTCLPGPCASTGPAAPLGAACGVGGPGAGQAASSQRGVCVLLPQESEAASRRSSPGWRRRP
NCGIRLPRIIRWRSVRQKRTQIPETRKRKDKGKGREPWSPTLWPPGLLG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 20.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_076420](#)

Locus ID: 65988


UniProt ID: [Q9BV97](#), [A0A024QZG5](#)

RefSeq Size: 2392



[View online »](#)

Cytogenetics: 16p11.2
RefSeq ORF: 573
Protein Families: Transcription Factors

Product images:

Coomassie blue staining of purified ZNF747 protein (Cat# TP300780). The protein was produced from HEK293T cells transfected with ZNF747 cDNA clone (Cat# [RC200780]) using MegaTran 2.0 (Cat# [TT210002]).