

Product datasheet for **TP300784**

NDUFB7 (NM_004146) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 7, 18kDa (NDUFB7), nuclear gene encoding mitochondrial protein, 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC200784 protein sequence
Red=Cloning site **Green**=Tags(s)

MG AHLVRRYLGDASVEPDPLQMPTFPPDYGFPERKEREMVATQQEMMDAQLRLQLRDYCAHHLIRLLKCK
RDSFPNFLACKQERHDWDYCEHRDYVMRMKEFERERLLQRKKRREKAAELAKGQGPGEVDPKVAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 16.2 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_004137](#)

Locus ID: 4713

UniProt ID: [P17568](#)

RefSeq Size: 573



[View online »](#)

Cytogenetics: 19p13.12


RefSeq ORF: 411

Synonyms: B18; CI-B18

Summary: The protein encoded by this gene is a subunit of the multisubunit NADH:ubiquinone oxidoreductase (complex I). Mammalian complex I is composed of 45 different subunits. It is located at the mitochondrial inner membrane. This protein has NADH dehydrogenase activity and oxidoreductase activity. It transfers electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone. [provided by RefSeq, Jul 2008]

Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

Product images:



Coomassie blue staining of purified NDUFB7 protein (Cat# TP300784). The protein was produced from HEK293T cells transfected with NDUFB7 cDNA clone (Cat# [RC200784]) using MegaTran 2.0 (Cat# [TT210002]).