

## Product datasheet for **TP300791L**

### **EFCAB11 (NM\_145231) Human Recombinant Protein**


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 14 open reading frame 143 (C14orf143), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200791 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MFFSEARARSRTWEASPSEHRKWVEVFKACDEDHKGYSREDFKTAVVMLFGYKPSKIEVDSVMSSINPN TSGILLEGFLNIVRKKKEAQRYRNEVRHIFTAFDTYYRGFLTLEDFKKAFRQVAPKLPERTVLEVFREVD RSDSGHVSFRDFEYALNYGQKEA
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	19 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_660274</a>
Locus ID:	90141
UniProt ID:	<a href="#">Q9BUY7</a>
RefSeq Size:	2868



[View online »](#)

Cytogenetics: 14q32.11  
RefSeq ORF: 489  
Synonyms: C14orf143

**Product images:**

Coomassie blue staining of purified EFCAB11 protein (Cat# [TP300791]). The protein was produced from HEK293T cells transfected with EFCAB11 cDNA clone (Cat# [RC200791]) using MegaTran 2.0 (Cat# [TT210002]).