

## Product datasheet for **TP300796M**

### MANBAL (NM\_001003897) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human mannosidase, beta A, lysosomal-like (MANBAL), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200796 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MASDLDFSPPEVPEPTFLENLLRYGLFLGAIQQLICVLAIIVPIPKSHEAEAEPSSEPRSAEVTRKPKAAV  
PSVNRKPKKETKKKR

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	9.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001003897</a>
Locus ID:	63905
UniProt ID:	<a href="#">Q9NQG1</a>
RefSeq Size:	1268




[View online »](#)

Cytogenetics: 20q11.23

RefSeq ORF: 255

Protein Families: Transmembrane

### Product images:



Coomassie blue staining of purified MANBAL protein (Cat# [TP300796]). The protein was produced from HEK293T cells transfected with MANBAL cDNA clone (Cat# [RC200796]) using MegaTran 2.0 (Cat# [TT210002]).