

Product datasheet for **TP300811M**

FSTL3 (NM_005860) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human follistatin-like 3 (secreted glycoprotein) (FSTL3), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200811 protein sequence Red =Cloning site Green =Tags(s)
	<p>MRPGAPGPLWPLPWGALAWAVGVSSMSGNPAPGGVCWLQQGQEATCSLVLQTDVTRAECASGNIDTA WSNLTHPGNKINLLGFLGLVHCLPCKDSCDGVCEGPGKACRMLGGRPRCECAPDCSGLPARLQVCGSDGA TYRDECELRAARCRGHPDLSVMYRGRCKRKSCEHVVCPRPQSCVVDQTGSAHCVVCRAAPCPVPSSPGQEL CGNNNVTYISSCHMRQATCFLGRSIGVRHAGSCAGTPEEPPGGESAEEEENFV</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	24.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005851
Locus ID:	10272
UniProt ID:	O95633 , A0A024R1Y8



[View online >](#)

RefSeq Size:	2525
Cytogenetics:	19p13.3
RefSeq ORF:	789
Synonyms:	FLRG; FSRP
Summary:	Follistatin-like 3 is a secreted glycoprotein of the follistatin-module-protein family. It may have a role in leukemogenesis. [provided by RefSeq, Jul 2008]
Protein Families:	Secreted Protein

Product images:



Coomassie blue staining of purified FSTL3 protein (Cat# [TP300811]). The protein was produced from HEK293T cells transfected with FSTL3 cDNA clone (Cat# [RC200811]) using MegaTran 2.0 (Cat# [TT210002]).