

Product datasheet for **TP300826**

DDA1 (NM_024050) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human DET1 and DDB1 associated 1 (DDA1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200826 protein sequence Red =Cloning site Green =Tags(s)
	MADFLKGLPVYNKSNFSRFHADSVCKASNRRPSVYLP TREY PSEQIIVTEKTNILLRYLHQW DKKNAAK KRDQEQVELEGESSAPPRKVARTDSPDMHEDT
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	11.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_076955
Locus ID:	79016
UniProt ID:	Q9BW61 , A0A024R7J7
RefSeq Size:	4078
Cytogenetics:	19p13.11




[View online »](#)

RefSeq ORF: 306

Synonyms: C19orf58; PCIA1

Summary: May be involved in ubiquitination and subsequent proteasomal degradation of target proteins. Component of the DDD-E2 complexes which may provide a platform for interaction with CUL4A and WD repeat proteins.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified DDA1 protein (Cat# TP300826). The protein was produced from HEK293T cells transfected with DDA1 cDNA clone (Cat# [RC200826]) using MegaTran 2.0 (Cat# [TT210002]).