

Product datasheet for **TP300879L**

PUS1 (NM_025215) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human pseudouridylate synthase 1 (PUS1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200879 protein sequence Red =Cloning site Green =Tags(s) MGLQLRALLGAFGRWTLRLGPRPSCSPRMAGNAEPPPAGAACPQDRRSCSGRAGGDRVWEDGEHPAKKLLK SGGDEERREKPPKRKIVLLMAYSGKGYHGMQRNVGSSQFKTIEDDLVSALVRSGCIPENHGEDMRKMSFQ RCARTDKGVSAAGQVWSLKVWLIDDILEKINSHLPSHIRILGLKRVTTGGFNSKNRCDARTYCYLLPTFAF AHKDRDVQDETYRLSAETLQQVNRLLACYKGTHNFHNFTSQKGPQDPSACRYILEMYCEEFPVREGLEFA VIRVKGQSFMMHQIRKMGVGLVAIVKGYAPESVLERSWGTEKVDVVPKAPGLGLVLERVHFEKYNQRFQND GLHEPLDWAQEEGKVAAFKEEHYPTIIGTERDERSMAQWLSTLPIHNFSATALTAGGTGAKVPSPLEGS EGDGD TD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	47.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_079491</u>



[View online »](#)

Locus ID:	80324
UniProt ID:	Q9Y606 , E5KMT5
RefSeq Size:	2014
Cytogenetics:	12q24.33
RefSeq ORF:	1281
Synonyms:	MLASA1
Summary:	This gene encodes a pseudouridine synthase that converts uridine to pseudouridine once it has been incorporated into an RNA molecule. The encoded enzyme may play an essential role in tRNA function and in stabilizing the secondary and tertiary structure of many RNAs. A mutation in this gene has been linked to mitochondrial myopathy and sideroblastic anemia. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Sep 2009]

Product images:



Coomassie blue staining of purified PUS1 protein (Cat# [TP300879]). The protein was produced from HEK293T cells transfected with PUS1 cDNA clone (Cat# [RC200879]) using MegaTran 2.0 (Cat# [TT210002]).