

Product datasheet for **TP300882M**

SARG (C1orf116) (NM_023938) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 1 open reading frame 116 (C1orf116), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200882 representing NM_023938 Red =Cloning site Green =Tags(s)
	MPERELWPAGTGSEPVTRVGSCDSMMSSTSTRSGSSDSSYDFLSTEEKECLLFLEETIGSLDTEADSGLS TDESEPATTPRGFRALPITQPTPRGGPEETITQQGRTPRTVTESSSSHPPEPQGLGLRSGSYSLPRNIHI ARSQNFRKSTTQASSHNPGEPEGRLAPEPEKEQVSQSSQPRQAPASPQEAALDLVDLIPPPEAFRDTQPE QCREASLPEGPGQQGHTPQLHTPSSSQEREQTPSEAMSQKAKETVSTRYTQPQPPPAGLPQNARAEDAP L SSGEDPNSRLAPLTTPKPRKLPPNIVLKSSRSSFHSDPQHWLSRHTEAAPGDSGLISCSLQEQRKARKEA LEKLGLPQDQDEPGLHLSKPTSSIRPKETRAQHLSAPAGLAQPAAPAQASAAIPAAAGKALAQAAPAPAGP AQGPLPMKSPAPGNVAASKSMPIPIKAPRANSALTTPKPESGLTLQESNTPGLRQMNFKSNTLERSGVG LSSYLSTEKDASPKTSTSLGKGSFLDKISPSVLNRNSRPRASLTGKDFAGIQVGKLADLEQEQSSKRLS YQGQSRDKLPRPPCVSVKISPKGVPNEHRREALKKLGLLKE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	63.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_076427</u>
Locus ID:	79098
UniProt ID:	<u>Q9BW04</u>
RefSeq Size:	5490
Cytogenetics:	1q32.1
RefSeq ORF:	1803
Synonyms:	SARG
Summary:	Putative androgen-specific receptor.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified C1orf116 protein (Cat# [TP300882]). The protein was produced from HEK293T cells transfected with C1orf116 cDNA clone (Cat# [RC200882]) using MegaTran 2.0 (Cat# [TT210002]).