

Product datasheet for **TP300933M**

RUVBL2 (NM_006666) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RuvB-like 2 (E. coli) (RUVBL2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200933 protein sequence Red =Cloning site Green =Tags(s)

MATVTATTKVPEIRDVTRIERIGAHSHIRGLGLDDALEPRQASQGMVGLAARRAAGVVLEMIREGKIAG
RAVLIAGQP GTGKTAIAMGMAQALGPDTPFTAAGSEIFSLEMSKTEALTQAFRRSIGVRIKEETEIIIEG
EVVEIQIDRPATGTGSKVGLTLKTTEMETIYDLGTKMIESLTKDKVQAGDVITIDKATGKISKLGRSFT
RARDYDAMGSQTKFVQCPDGELQKRKEVHTVSLHEIDVINSRTQGFALFSGDTGEIKSEVREQINAKV
AEWREEGKAEIIPGVLFIDEVHMLDIESFSFLNRALES DMAPVLIMATNRGITRIGTSYQSPHGIPIDL
LDRLIVSTTPYSEKDTKQILRIRCEEEEDVEMSEDAYTVLTRIGLETSLRYAIQLITAASLVCRKRKGTGTE
VQVDDIKRVYSLFLDESRSTQYMKEYQDAFLFNLKGETMDTS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	51 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_006657</u>



[View online »](#)

Locus ID:	10856
UniProt ID:	Q9Y230
RefSeq Size:	1488
Cytogenetics:	19q13.33
RefSeq ORF:	1389
Synonyms:	CGI-46; ECP-51; ECP51; INO80J; REPTIN; RVB2; TAP54-beta; TIH2; TIP48; TIP49B
Summary:	This gene encodes the second human homologue of the bacterial RuvB gene. Bacterial RuvB protein is a DNA helicase essential for homologous recombination and DNA double-strand break repair. Functional analysis showed that this gene product has both ATPase and DNA helicase activities. This gene is physically linked to the CGB/LHB gene cluster on chromosome 19q13.3, and is very close (55 nt) to the LHB gene, in the opposite orientation. [provided by RefSeq, Jul 2008]
Protein Families:	Transcription Factors

Product images:



Coomassie blue staining of purified RUVBL2 protein (Cat# [TP300933]). The protein was produced from HEK293T cells transfected with RUVBL2 cDNA clone (Cat# [RC200933]) using MegaTran 2.0 (Cat# [TT210002]).