

Product datasheet for **TP300944**

PAXX (NM_183241) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 9 open reading frame 142 (C9orf142), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200944 protein sequence Red =Cloning site Green =Tags(s)
	MDPLSPPLCTLPPGPEPPRFVFCYCEGEESGEGDRGGFNLYVTDAELWSTCFTPDSLAAALKARFGLSAAE DITPRFRAACEQQAVALTQEDRASLTLSGGPSALAFDLSKVPGEAAPRLRALTLGLAKRVWSLERRLA AAEETAVSPRKSPRPAGPQLFLPDPDPQRRGGPGVRRRCPGESLINPGFKSKKPAGGVDFDET
	SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_899064
Locus ID:	286257
UniProt ID:	Q9BUH6
RefSeq Size:	851



[View online »](#)


Cytogenetics: 9q34.3

RefSeq ORF: 612

Synonyms: C9orf142; XLS

Summary: The protein encoded by this gene plays a role in the nonhomologous end joining (NHEJ) pathway of DNA double-strand break repair. The encoded protein may function to stabilize the Ku70/Ku80 heterodimer to facilitate the assembly and maintain the stability of the NHEJ complex. [provided by RefSeq, Jul 2016]

Product images:



Coomassie blue staining of purified PAXX protein (Cat# TP300944). The protein was produced from HEK293T cells transfected with PAXX cDNA clone (Cat# [RC200944]) using MegaTran 2.0 (Cat# [TT210002]).