

Product datasheet for **TP300951L**

LRR61 (NM_023942) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human leucine rich repeat containing 61 (LRR61), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200951 protein sequence Red =Cloning site Green =Tags(s)

MDPPAEKPGEAGGLQITPQLLKSRTGEFSLESILLKLRGLGLADLGLGCEGLGLEWLDLSGNALHGLP
LASLRQLAVLNVSNNRLTGLEPLATCENLQSLNAAGNLLATPGQLQCLAGLPCLEYLRLRDPLARLSNPL
CANPSYWAAVRELLPGLKVIDGERVIGRGSEFYQLCRDLSSLRPSSSPGPRATEAQWVEPGYWESWPS
RSSSILEEACRQFQDTLQECWDLDRQASDSLAQAEQVLSSAGPTSSVF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	27.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_076431</u>
Locus ID:	65999



[View online »](#)

UniProt ID:	<u>Q9BV99</u>
RefSeq Size:	1809
Cytogenetics:	7q36.1
RefSeq ORF:	777
Synonyms:	HSPC295

Product images:

Coomassie blue staining of purified LRRC61 protein (Cat# [TP300951]). The protein was produced from HEK293T cells transfected with LRRC61 cDNA clone (Cat# [RC200951]) using MegaTran 2.0 (Cat# [TT210002]).