

Product datasheet for **TP300955**

MXRA8 (NM_032348) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human matrix-remodelling associated 8 (MXRA8), 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC200955 protein sequence
Red=Cloning site **Green**=Tags(s)

MALPSRILLWKLVLQSSAVLLHSGSSVPAAGSSVSESAVSWEAGARAVLRQCSPRMVWTQDRLHDR
Q
RVLHWDLRGPGGGPARRLLDLYSAGEQRVYEARDRGRLELSASAFDDGNFSLIRAVEETDAGLYTCNLH
HHYCHLYESLAVRLEVDGPPATPAYWDGEKEVLAVARGAPALLTCVNRGHVWTDHRVVEEAQQVHWD
RQ
PPGVPHDRADRLLDLYASGERRAYGPLFLRDRVAVGADAFERGFSLRIEPLVADEGTYSCHLHHHYCG
LHERRVFHLTVAEPHAEPPIRGSPGNGSSHSGAPGPDPTLARGHNVINVIVPESRAHFFQQLGYVLATLL
LFILLVTVLLAARRRRGGYEYSDQKSGKSKGKDVNLAEFAVAAGDQMLYRSEDIQLDYKNILKERAEL
AHSPLPAKYIDLKDGFRKENCK

LEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 49 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


Storage: Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_115724
Locus ID:	54587
UniProt ID:	Q9BRK3
RefSeq Size:	2534
Cytogenetics:	1p36.33
RefSeq ORF:	1326
Synonyms:	ASP3
Summary:	<p>Transmembrane protein which can modulate activity of various signaling pathways, probably via binding to integrin ITGAV:ITGB3 (PubMed:22492581, PubMed:23386276). Mediates heterophilic cell-cell interactions in vitro (By similarity). Inhibits osteoclastogenesis downstream of TNFSF11/RANKL and CSF1, where it may function by attenuating signaling via integrin ITGB3 and MAP kinase p38 (By similarity). Plays a role in cartilage formation where it promotes proliferation and maturation of growth plate chondrocytes (By similarity). Stimulates formation of primary cilia in chondrocytes (By similarity). Enhances expression of genes involved in the hedgehog signaling pathway in chondrocytes, including the hedgehog signaling molecule IHH; may also promote signaling via the PTHLH/PTHrP pathway (By similarity). Plays a role in angiogenesis where it suppresses migration of endothelial cells and also promotes their apoptosis (PubMed:23386276). Inhibits VEGF-induced activation of AKT and p38 MAP kinase in endothelial cells (PubMed:23386276). Also inhibits VTN (vitronectin)-mediated integrin ITGAV:ITGB3 signaling and activation of PTK2/FAK (PubMed:23386276). May play a role in the maturation and maintenance of the blood-brain barrier (By similarity). [UniProtKB/Swiss-Prot Function]</p>
Protein Families:	Transmembrane

Product images:



Coomassie blue staining of purified MXRA8 protein (Cat# TP300955). The protein was produced from HEK293T cells transfected with MXRA8 cDNA clone (Cat# [RC200955]) using MegaTran 2.0 (Cat# [TT210002]).