

Product datasheet for **TP300957**

Cytochrome C Oxidase subunit VIb (COX6B1) (NM_001863) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cytochrome c oxidase subunit VIb polypeptide 1 (ubiquitous) (COX6B1), nuclear gene encoding mitochondrial protein
Species:	Human
Expression Host:	HEK293T
Tag:	C-Myc/DDK
Predicted MW:	10 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001854
Locus ID:	1340
RefSeq Size:	590
Cytogenetics:	19q13.12
RefSeq ORF:	258
Synonyms:	COX6B; COXG; COXVIb1




[View online »](#)

Summary:

Cytochrome c oxidase (COX), the terminal enzyme of the mitochondrial respiratory chain, catalyzes the electron transfer from reduced cytochrome c to oxygen. It is a heteromeric complex consisting of 3 catalytic subunits encoded by mitochondrial genes and multiple structural subunits encoded by nuclear genes. The mitochondrially-encoded subunits function in electron transfer, and the nuclear-encoded subunits may be involved in the regulation and assembly of the complex. This nuclear gene encodes subunit VIb. Mutations in this gene are associated with severe infantile encephalomyopathy. Three pseudogenes COX6BP-1, COX6BP-2 and COX6BP-3 have been found on chromosomes 7, 17 and 22q13.1-13.2, respectively. [provided by RefSeq, Jan 2010]

Protein Pathways:

Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

Product images:

Coomassie blue staining of purified COX6B1 protein (Cat# TP300957). The protein was produced from HEK293T cells transfected with COX6B1 cDNA clone (Cat# [RC200957]) using MegaTran 2.0 (Cat# [TT210002]).