

Product datasheet for **TP300962**

GGCT (NM_024051) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human gamma-glutamyl cyclotransferase (GGCT), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200962 protein sequence Red =Cloning site Green =Tags(s) MANS GCKDVTGPDEESFLYFAYGSNLLTERIHLRNPSAAFFCVARLQDFKLDFGNSQGKTSQ TWHGGIAT IFQSPGDEVWGVWKMNKS NLSLDEQEGVKSGMYV VIEVKVATQEGKEITCRSYLMTNYESAPPSPQYK KIICMGAKENGLPLEYQEK LKAIEPN DYTGKVSEEIEDIIKKGETQTL TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_076956</u>
Locus ID:	79017
UniProt ID:	<u>O75223</u>
RefSeq Size:	1197


[View online »](#)

Cytogenetics:	7p14.3
RefSeq ORF:	564
Synonyms:	C7orf24; CRF21; GCTG; GGC
Summary:	The protein encoded by this gene catalyzes the formation of 5-oxoproline from gamma-glutamyl dipeptides, the penultimate step in glutathione catabolism, and may play a critical role in glutathione homeostasis. The encoded protein may also play a role in cell proliferation, and the expression of this gene is a potential marker for cancer. Pseudogenes of this gene are located on the long arm of chromosome 5 and the short arm of chromosomes 2 and 20. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2010]

Protein Pathways: Glutathione metabolism

Product images:



Coomassie blue staining of purified GGCT protein (Cat# TP300962). The protein was produced from HEK293T cells transfected with GGCT cDNA clone (Cat# [RC200962]) using MegaTran 2.0 (Cat# [TT210002]).