

Product datasheet for **TP300977L**

SF3B5 (NM_031287) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human splicing factor 3b, subunit 5, 10kDa (SF3B5), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200977 protein sequence Red =Cloning site Green =Tags(s)
	MTDRYTIHSQLEHLQSKYIGTGHADTTKWEVLVNQHRDSYCSYMGHFDLLNYFAIAENESKARVRFNLME KMLQPCGPPADKPEEN
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	10 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_112577
Locus ID:	83443
UniProt ID:	Q9BWJ5
RefSeq Size:	786
Cytogenetics:	6q24.2



[View online »](#)


RefSeq ORF: 258

Synonyms: SF3b10; Ysf3

Summary: Involved in pre-mRNA splicing as a component of the splicing factor SF3B complex, a constituent of the spliceosome (PubMed:27720643, PubMed:28781166). SF3B complex is required for 'A' complex assembly formed by the stable binding of U2 snRNP to the branchpoint sequence (BPS) in pre-mRNA. Sequence independent binding of SF3A/SF3B complex upstream of the branch site is essential, it may anchor U2 snRNP to the pre-mRNA (PubMed:12234937).[UniProtKB/Swiss-Prot Function]

Protein Pathways: Spliceosome

Product images:



Coomassie blue staining of purified SF3B5 protein (Cat# [TP300977]). The protein was produced from HEK293T cells transfected with SF3B5 cDNA clone (Cat# [RC200977]) using MegaTran 2.0 (Cat# [TT210002]).