

Product datasheet for **TP300979**

JAB1 (COPS5) (NM_006837) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human COP9 constitutive photomorphogenic homolog subunit 5 (Arabidopsis) (COPS5), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200979 protein sequence Red =Cloning site Green =Tags(s)

MAASGSGMAQKTWELANNMQEAQSIDEIYKYDKKQQQEILAAKPWTKDHHYFKYCKISALALLKMVMHAR
SGGNLEVMGLMLGKVDGETMIIMDSFALPVEGTETRVNAQAAAYEYMAAYIENAKQVGRLENAIGWYHSH
PGYGCWLSGIDVSTQMLNQQFQEPFVAWIDPTRTISAGKVNLGAFRTYPKGYKPPDEGPSEYQTIPLNK
IEDFGVHCKQYYALEVSYFKSSLDRKLELLWNKYWVNTLSSSSLLTNADYTTGQVFDLSEKLEQSEAQL
GRGSFMLGLETHDRKSEDKLAKATRDSCKTIEAIHGLMSQVIKDKLFNQINIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	37.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006828
Locus ID:	10987



[View online »](#)

UniProt ID: [Q92905](#), [A0A024R7W9](#)

RefSeq Size: 1510

Cytogenetics: 8q13.1


RefSeq ORF: 1002

Synonyms: CSN5; JAB1; MOV-34; SGN5

Summary: The protein encoded by this gene is one of the eight subunits of COP9 signalosome, a highly conserved protein complex that functions as an important regulator in multiple signaling pathways. The structure and function of COP9 signalosome is similar to that of the 19S regulatory particle of 26S proteasome. COP9 signalosome has been shown to interact with SCF-type E3 ubiquitin ligases and act as a positive regulator of E3 ubiquitin ligases. This protein is reported to be involved in the degradation of cyclin-dependent kinase inhibitor CDKN1B/p27Kip1. It is also known to be an coactivator that increases the specificity of JUN/AP1 transcription factors. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Protease, Transcription Factors

Product images:



Coomassie blue staining of purified COPS5 protein (Cat# TP300979). The protein was produced from HEK293T cells transfected with COPS5 cDNA clone (Cat# [RC200979]) using MegaTran 2.0 (Cat# [TT210002]).