

Product datasheet for **TP300988**

SSBP4 (NM_032627) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human single stranded DNA binding protein 4 (SSBP4), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC200988 protein sequence Red =Cloning site Green =Tags(s) MYAKGGKGSAPVSDSQAREKLALYVYELLHIGAQKSAQTFLEIRWEKNITLGEPPGFLHSWWCVFWDLYCAAPDRREACEHSGEAKAFQDYSAAAAPSPVMGSMAPGDTMAAGSMAAGFFQGPPGSQPSHPNPNAPMMGPHGQPFMSPRFPGGPRPTLRMPSPQPPAGLPGSQPLLPGAMEPSRAQGHPSMGGPMQRVTPPRGMASVGPQSYGGMRPPNSLAGPGLPAMNMGPGVRGPWASPSGNSIPYSSSSPGSYTGPPGGGGPPGTPIMPSPGDSTNSSENMYTIMNPIGQAGRANFPLGPGPEGPMAAMSAMEPHHVNGSLGSGDMDGLPKSSPGAVAGLSNAPGTPRDDGEMAAAGTFLHPFPSESYSPGMTMSV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	39.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_116016</u>



[View online »](#)

Locus ID: 170463
UniProt ID: [Q9BWG4](#)
RefSeq Size: 1817
Cytogenetics: 19p13.11
RefSeq ORF: 1155
Protein Families: Transcription Factors

Product images:



Coomassie blue staining of purified SSBP4 protein (Cat# TP300988). The protein was produced from HEK293T cells transfected with SSBP4 cDNA clone (Cat# [RC200988]) using MegaTran 2.0 (Cat# [TT210002]).