

Product datasheet for **TP301032M**

TRIAP1 (NM_016399) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human TP53 regulated inhibitor of apoptosis 1 (TRIAP1), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201032 protein sequence Red =Cloning site Green =Tags(s)
	MNSVGEACTDMKREYDQCFNRWFAEKFLKGDSSGDPCTDLFKRYQQCVQKAIKEKEIPIEGLEFMGHGKE KPENSS
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	8.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057483
Locus ID:	51499
UniProt ID:	Q43715
RefSeq Size:	1196
Cytogenetics:	12q24.31




[View online »](#)

RefSeq ORF: 228

Synonyms: HSPC132; MDM35; P53CSV; WF-1

Summary: Involved in the modulation of the mitochondrial apoptotic pathway by ensuring the accumulation of cardiolipin (CL) in mitochondrial membranes. In vitro, the TRIAP1:PRELID1 complex mediates the transfer of phosphatidic acid (PA) between liposomes and probably functions as a PA transporter across the mitochondrion intermembrane space to provide PA for CL synthesis in the inner membrane (PubMed:23931759). Likewise, the TRIAP1:PRELID3A complex mediates the transfer of phosphatidic acid (PA) between liposomes (in vitro) and probably functions as a PA transporter across the mitochondrion intermembrane space (in vivo) (PubMed:26071602). Mediates cell survival by inhibiting activation of caspase-9 which prevents induction of apoptosis (PubMed:15735003).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified TRIAP1 protein (Cat# [TP301032]). The protein was produced from HEK293T cells transfected with TRIAP1 cDNA clone (Cat# [RC201032]) using MegaTran 2.0 (Cat# [TT210002]).