

Product datasheet for **TP301047M**

CENPM (NM_024053) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human centromere protein M (CENPM), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201047 protein sequence Red =Cloning site Green =Tags(s)
	 MSVLRPLDKLPGLNTATILLVGTEDALLQQLADSMKEDCASELKVHLAKSLPLPSSVNRPRIDLIVFVV NLHISKYSLQNTTEESLRHVDASFFLGKVCFLATGAGRESHCSHRHTVVKLAHTYQSPLLYCDLEVEGFRA TMAQRLVRVLQICAGHVPGVSA LNLLSLLRSSEGPSLEDL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	19.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_076958
Locus ID:	79019
UniProt ID:	Q9NSP4 , A0A024R1Q3
RefSeq Size:	947



[View online »](#)

Cytogenetics: 22q13.2


RefSeq ORF: 540

Synonyms: C22orf18; CENP-M; PANE1

Summary: The protein encoded by this gene is an inner protein of the kinetochore, the multi-protein complex that binds spindle microtubules to regulate chromosome segregation during cell division. It belongs to the constitutive centromere-associated network protein group, whose members interact with outer kinetochore proteins and help to maintain centromere identity at each cell division cycle. The protein is structurally related to GTPases but cannot bind guanosine triphosphate. A point mutation that affects interaction with another constitutive centromere-associated network protein, CENP-I, impairs kinetochore assembly and chromosome alignment, suggesting that it is required for kinetochore formation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2015]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified CENPM protein (Cat# [TP301047]). The protein was produced from HEK293T cells transfected with CENPM cDNA clone (Cat# [RC201047]) using MegaTran 2.0 (Cat# [TT210002]).