

Product datasheet for TP301077

IVD (NM_002225) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human isovaleryl Coenzyme A dehydrogenase (IVD), nuclear gene encoding mitochondrial protein, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC201077 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MATATRLLGCRVASWRLRPPLAGFVSQRAHSLLPVDDAINGLSEEQRQLRQTMAKFLQEHLAPKAQEIDR
SNEFKNLREFWKQLGNLGVLGITAPVQYGGSLGYLEHVLVMEEISRASGAVGLSYGAHSNLCINQLVRN
GNEAQKEKYLPLKISGEYIGALAMSEPNAAGSDVSMKCLKAEKGNHYILNGNKFWITNGPDADVLIYAK
TDLAAVPASRGITAFIVEKGMPPGFSTSKKLDKLGMRGNTCELIFEDCKIPAANILGHENKGVYVLMVSL
DLERLVLGGPLGLMQAVLDHTIPLYLVREAFGQKIGHFQLMQGKMADMYTRLMACRQYVYNVAKACDEG
HCTAKDCAGVILYSAECATQVALDGIQCFGGNGYINDFPMGRFLRDAKLYEIGAGTSEVRRLLVIGRAFNA
DFH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	46.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_002216](#)

Locus ID: 3712

UniProt ID: [P26440](#), [A0A0A0MT83](#)

RefSeq Size: 4673

Cytogenetics: 15q15.1

RefSeq ORF: 1269


Synonyms: ACAD2; IVDH

Summary: Isovaleryl-CoA dehydrogenase (IVD) is a mitochondrial matrix enzyme that catalyzes the third step in leucine catabolism. The genetic deficiency of IVD results in an accumulation of isovaleric acid, which is toxic to the central nervous system and leads to isovaleric acidemia. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Aug 2017]

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Valine, leucine and isoleucine degradation

Product images:



Coomassie blue staining of purified IVD protein (Cat# TP301077). The protein was produced from HEK293T cells transfected with IVD cDNA clone (Cat# [RC201077]) using MegaTran 2.0 (Cat# [TT210002]).