

Product datasheet for TP301085L

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

NSE (ENO2) (NM_001975) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human enolase 2 (gamma, neuronal) (ENO2), 1 mg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC201085 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s)

MSIEKIWAREILDSRGNPTVEVDLYTAKGLFRAAVPSGASTGIYEALELRDGDKQRYLGKGVLKAVDHIN STIAPALISSGLSVVEQEKLDNLMLELDGTENKSKFGANAILGVSLAVCKAGAAERELPLYRHIAQLAGN SDLILPVPAFNVINGGSHAGNKLAMQEFMILPVGAESFRDAMRLGAEVYHTLKGVIKDKYGKDATNVGDE GGFAPNILENSEALELVKEAIDKAGYTEKIVIGMDVAASEFYRDGKYDLDFKSPTDPSRYITGDQLGALY QDFVRDYPVVSIEDPFDQDDWAAWSKFTANVGIQIVGDDLTVTNPKRIERAVEEKACNCLLLKVNQIGSV TEAIQACKLAQENGWGVMVSHRSGETEDTFIADLVVGLCTGQIKTGAPCRSERLAKYNQLMRIEEELGDE

ARFAGHNFRNPSVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 47.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by

conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001966



Locus ID: 2026

UniProt ID: <u>P09104</u>, <u>Q6FHV6</u>

RefSeq Size: 2423

Cytogenetics: 12p13.31

RefSeq ORF: 1302

Synonyms: HEL-S-279; NSE


Summary: This gene encodes one of the three enclase isoenzymes found in mammals. This isoenzyme, a

homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

[provided by RefSeq, Jul 2008]

Protein Pathways: Glycolysis / Gluconeogenesis, Metabolic pathways, RNA degradation

Product images:



Coomassie blue staining of purified ENO2 protein (Cat# [TP301085]). The protein was produced from HEK293T cells transfected with ENO2 cDNA clone (Cat# [RC201085]) using MegaTran 2.0 (Cat# [TT210002]).