

Product datasheet for TP301092L

SERPBP1 (NM_001018067) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SERPINE1 mRNA binding protein 1 (SERBP1), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201092 representing NM_001018067 Red =Cloning site Green =Tags(s)
	<p>MPGHLQEGFGCVTNRFDQLFDDSDPFEVLKAAENKKKEAGGGGVGGPAGKSAQAQAQTNSNAAGKQL RKESQKDRKNPLPPSVGVVDKKEETQPPVALKKEGIRRVGRRPDQQLQGEGKIIDRRPERRPPRERRFEK PLEEKGEGGFSVDRPIIDRPIRGRGGLGRGRGGRGRGMGRGDGFDSRGKREFDRHSGSDRSSFSHYSGL KHEDKRGGSGSHNWGTVKDELTEPKYIQKQISYNYSDLDQSNVTEETPEGEEHHPVADTENKENEVEEV KEEGPKEMTLDEWKAIQNKDRAKVEFNIRKPNEGADGQWKKGFVLHKSSEEHAEDSVMDHHRKFPAND ITSQLEINFGDLGRPRGGRGGRGGRGGRGGRPNRGRSRTDKSSASAPDVEDDPEAFPALA</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	44.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001018077</u>



[View online »](#)

Locus ID:	26135
UniProt ID:	Q8NC51
RefSeq Size:	6764
Cytogenetics:	1p31.3
RefSeq ORF:	1224
Synonyms:	CGI-55; CHD3IP; HABP4L; PAI-RBP1; PAIRBP1
Summary:	May play a role in the regulation of mRNA stability. Binds to the 3'-most 134 nt of the SERPINE1/PAI1 mRNA, a region which confers cyclic nucleotide regulation of message decay. Seems to play a role in PML-nuclear bodies formation (PubMed:28695742).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified SERBP1 protein (Cat# [TP301092]). The protein was produced from HEK293T cells transfected with SERBP1 cDNA clone (Cat# [RC201092]) using MegaTran 2.0 (Cat# [TT210002]).