

Product datasheet for **TP301103L**

ARPC1B (NM_005720) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human actin related protein 2/3 complex, subunit 1B, 41kDa (ARPC1B), 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC201103 protein sequence
Red=Cloning site **Green**=Tags(s)

MAYHSFLVEPISCHAWNKDRTQIAICPNNHEVHIYEKSGAKWTKVHELKEHNGQVTGIDWAPESNRIVTC
GTDNRNAYVWTLKGRWTKPTLVILRINRAARCVRWAPNENKFAVSGSRVISICYFEQENDWWVCKHIKPK
IRSTVLSLDWHPNPNVLLAAGSCDFKCRIFSAYIKEVEERPAPTWPWGSKMPFGELMFESSSSCGWVHGVCF
SASGSRVAWVSHDSTVCLADADKKMAVATLASETLPLLALTFITDNSLVAAGHDCFPVLFTYDAAAGMLS
FGGRLDVPKQSSQRGLTARERFQNLDDKKASSEGGTAAGAGLDSLHKNSVSVLSVLSGGKAKCSQFCTTGM
DGGMSIWDVKSLESALKDLKIK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 40.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_005711](#)



[View online »](#)

Locus ID: 10095

UniProt ID: [O15143](#), [A4D275](#)

RefSeq Size: 1551

Cytogenetics: 7q22.1


RefSeq ORF: 1116

Synonyms: ARC41; IMD71; p40-ARC; p41-ARC; PLTEID

Summary: This gene encodes one of seven subunits of the human Arp2/3 protein complex. This subunit is a member of the SOP2 family of proteins and is most similar to the protein encoded by gene ARPC1A. The similarity between these two proteins suggests that they both may function as p41 subunit of the human Arp2/3 complex that has been implicated in the control of actin polymerization in cells. It is possible that the p41 subunit is involved in assembling and maintaining the structure of the Arp2/3 complex. Multiple versions of the p41 subunit may adapt the functions of the complex to different cell types or developmental stages. This protein also has a role in centrosomal homeostasis by being an activator and substrate of the Aurora A kinase. [provided by RefSeq, Mar 2011]

Protein Pathways: Fc gamma R-mediated phagocytosis, Pathogenic Escherichia coli infection, Regulation of actin cytoskeleton

Product images:



Coomassie blue staining of purified ARPC1B protein (Cat# [TP301103]). The protein was produced from HEK293T cells transfected with ARPC1B cDNA clone (Cat# [RC201103]) using MegaTran 2.0 (Cat# [TT210002]).