

## Product datasheet for TP301107

### CCDC9 (NM\_015603) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 9 (CCDC9), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201107 representing NM_015603 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MAATLDLKSKEEKDAELDKRIEALRRKNEALIRRYQEIEEDRKKAIELEGVAVTAPRKGRSVEKENVAVES EKNLGPSSRRSPGTPRPPGASKGGRTPPQQGGRAGMGRASRSWEGSPGEQPRGGGAGGRRRGRGRGSPHL SGAGDTSISDRKSKEWEERRRQNIKMNEMEKIAEYERNQREGVLEPNPVRNFLDDPRRRSGPLEESER DRREESRRHGRNWGGPDFERVRCGLEHERQGRAGLGSAGDMTMSMTGRERSEYLRWKQEREKIDQERLQ RHRKPTGQWRREWDAEKTDMFKDGPVPAHEPSHRYDDQAWARPPKPTTFGEFLSQHKAESSRRRRKSS RPQAKAAPRAYSDHDDRWETKEGAASPAPETPQPTSPETSPKETPMQPPEIPAPAHPPPEDEGEENEGEE DEEWEDISEDEEEEEIEVEEGDEEPEAQDHQAPEAAPTGIPCSEQAHGVPFSPPEEPLLEPQAPGTPSSPF SPPSGHQPVSDWGEEVELNSPRTTHLAGALSPGEAWPFESV</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	59.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP\\_056418](#)

Locus ID: 26093


UniProt ID: [Q9Y3X0](#), [B4DXW2](#)

RefSeq Size: 2026

Cytogenetics: 19q13.32

RefSeq ORF: 1593

### Product images:



Coomassie blue staining of purified CCDC9 protein (Cat# TP301107). The protein was produced from HEK293T cells transfected with CCDC9 cDNA clone (Cat# [RC201107]) using MegaTran 2.0 (Cat# [TT210002]).