

Product datasheet for TP301110L

GLI1 (NM_005269) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human GLI family zinc finger 1 (GLI1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201110 protein sequence Red =Cloning site Green =Tags(s)

MFNSMTPPPISSYGEPCCLRPLPSQGAPSVGTEGLSGPPFCHQANLMSGPHSYGPARETNSCTEGPLFSS
PRSAVKLTKRALSISPLSDASLDLQTVIRTSPSSLVAFINSRCTSPGGSYGHLSIGTMSPSLGFPAQMN
HQKGPSPFSGVQPCGPHDSARGGMIPHPQSRGPFPTCQLKSELDMLVGKCREEPLEGDMSSPNSTGIQDP
LLGMLDGREDLEREEKREPESVYETDCRWDGCSQEFDSEQLVHHINSEHIHGERKEFVCHWGGCSREL
PFKAQYMLVVHMRRHTGEKPHKCTFEGCRKSYSRLENLKTHLRSHTGEKPYMCEHEGCSKAFSNASDRAK
HQNRTHSNEKPYVCKLPGCTKRYTDPSSLRKHVKTVHGPDAHVTKRHRGDGPLPRAPSISTVEPKREREG
GPIREESRLTVPEGAMKPQSPGAQSSCSDHSPAGSAANTDSGVEMTGNAGGSTEDLSSLDEGPCIAGT
GLSTLRLENLRLDQLHQLRPIGTRGLKPLSHTGTTVSRRVGPVSLERRSSSSSSISSAYTVSRRSS
LASPPPGSPENGASSLPGLMPAQHYLLRARYASARGGGTSPTAASSLDRIGGLMPPWRSRAEYPGYN
PNAGVTRRASDPAQAADRPAPARVQRFKSLGCVHTPPTVAGGGQNFDPYLPYLVSPQPPSITENAAMDA
RGLQEEPEVGTSMVGSGLNPYMDFPPTDTLGYGGPEGAAAEPYGARGPGLPLGPGPPTNYGPNPCPQQA
SYPDPTQETWGEFSPSHSGLYPGPKALGGTYSQCPRLHYGQVQVKPEQGCPVGSSTGLAPCLNAHPSEG
PPHPQPLFSHYPPQSPPPYQLQSGPYTQPPDYLPSEPRPCLDFDSTHSTGQLKAQLVCNYVQSQQELLW
EGGGREDAPAEPSYQSPKFLGGSQVSPSRAPVNTYGPFGPNLPHKSGSYPTSPCHENFVVGANR
ASHRAAAPRLLPPLPTCYGPLKVGGTNPSCGHPEVGRLGGGPALYPPPEGQVCNPLDSLDDNTQLDFV
AILDEPQGLSPPPSHDQRGSSGHTPPPSGPPNMAVGNMSVLLRSLPGETEFLNSSA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	117.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Bioactivity:	In vitro kinase assay substrate (PMID: 28735864) Binding assay (PMID: 30035333)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005260
Locus ID:	2735
UniProt ID:	P08151
RefSeq Size:	3618
Cytogenetics:	12q13.3
RefSeq ORF:	3318
Synonyms:	GLI; PAPA8; PPD1
Summary:	This gene encodes a member of the Kruppel family of zinc finger proteins. The encoded transcription factor is activated by the sonic hedgehog signal transduction cascade and regulates stem cell proliferation. The activity and nuclear localization of this protein is negatively regulated by p53 in an inhibitory loop. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009]
Protein Families:	Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Stem cell relevant signaling - DSL/Notch pathway, Transcription Factors
Protein Pathways:	Basal cell carcinoma, Hedgehog signaling pathway, Pathways in cancer

Product images:



Coomassie blue staining of purified GLI1 protein (Cat# [TP301110]). The protein was produced from HEK293T cells transfected with GLI1 cDNA clone (Cat# [RC201110]) using MegaTran 2.0 (Cat# [TT210002]).