

## Product datasheet for TP301144L

### Phosphoribosyl pyrophosphate amidotransferase (PPAT) (NM\_002703) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human phosphoribosyl pyrophosphate amidotransferase (PPAT), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201144 protein sequence Red=Cloning site Green=Tags(s)

MELEELGIREECGVFGCIASGEWPTQLDVPHVITLGLVGLQHRGQESAGIVTSDGSSVPTFKSHKGMGLV  
NHVFTEDNLKLYVSNLGIGHTRYATTGKCELENCQPFVETLHGKIAVAHNGELVNAARLRKLLRHGI  
GLTSSDSEMITQLLAYTPPQEQDDTPDWVARIKMLMKEAPTAYSLLIMHRDVIYAVRDPYGNRPLCIGR  
LIPVSDINDKEKKTSETEGWVVSSESCSFLSIGARYYREVLPGIVEISRHNVQTLDIISRSEGNPVAFC  
IFEYVVFARPD SMFEDQM VYTVRYRCGQQLAIEAPVDADLVSTVPESATPAALAYAGKCGLPYVEVLCKN  
RYVGRTFIQPNMRLRQLGVAKKFGVLSDNFKGKRIVLDDSIVRGNTISPIIKLLKESGAKEVHIRVASP  
PIKYP CFMGINIPTKEELIANKPEFDHLAEYLGANSVWYLSVEGLVSSVQEGIKFKKQKEKKHDIMI QEN  
GNGLECFEKS GHCTACTLGKYPVELEW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	57.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_002694</a>
<b>Locus ID:</b>	5471
<b>UniProt ID:</b>	<a href="#">Q06203</a> , <a href="#">A8K4H7</a> , <a href="#">Q59G63</a>
<b>RefSeq Size:</b>	3713
<b>Cytogenetics:</b>	4q12
<b>RefSeq ORF:</b>	1551
<b>Synonyms:</b>	ATASE; GPAT; PRAT
<b>Summary:</b>	The protein encoded by this gene is a member of the purine/pyrimidine phosphoribosyltransferase family. It is a regulatory allosteric enzyme that catalyzes the first step of de novo purine nucleotide biosynthetic pathway. This gene and PAICS/AIRC gene, a bifunctional enzyme catalyzing steps six and seven of this pathway, are located in close proximity on chromosome 4, and divergently transcribed from an intergenic region. [provided by RefSeq, Mar 2011]
<b>Protein Families:</b>	Druggable Genome, Protease
<b>Protein Pathways:</b>	Alanine, aspartate and glutamate metabolism, Metabolic pathways, Purine metabolism

### Product images:



Coomassie blue staining of purified PPAT protein (Cat# [TP301144]). The protein was produced from HEK293T cells transfected with PPAT cDNA clone (Cat# [RC201144]) using MegaTran 2.0 (Cat# [TT210002]).