

## Product datasheet for **TP301206M**

### **Sigma1 receptor (SIGMAR1) (NM\_005866) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human sigma non-opioid intracellular receptor 1 (SIGMAR1), transcript variant 1, 100 µg
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>RC201206 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MQWAVGRRWAWAALLLAVAAVLTQVWWLWLGTSFVQREEIAQLARQYAGLDHELAFSRLIVELRRRLHP GHVLPDEELQWVFNAGGWMGAMCLLHASLSEYVLLFGTALGSRGHSGRYWAEISDTIISGTFHQWREGT TKSEVFYPGETVHGPGEATAVEWGPNTWMVEYGRGVIPSTLAFALADTVFSTQDFLTLFYTLRSYARGL RLELTTYLFGQDP  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
<b>Tag:</b>	C-Myc/DDK
<b>Predicted MW:</b>	24.9 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<u><a href="#">NP_005857</a></u>
<b>Locus ID:</b>	10280



[View online »](#)

UniProt ID: [Q99720](#)

RefSeq Size: 1728

Cytogenetics: 9p13.3


RefSeq ORF: 669

Synonyms: ALS16; DSMA2; hSigmaR1; OPRS1; SIG-1R; sigma1R; SR-BP; SR-BP1; SRBP

**Summary:** This gene encodes a receptor protein that interacts with a variety of psychotomimetic drugs, including cocaine and amphetamines. The receptor is believed to play an important role in the cellular functions of various tissues associated with the endocrine, immune, and nervous systems. As indicated by its previous name, opioid receptor sigma 1 (OPRS1), the product of this gene was erroneously thought to function as an opioid receptor; it is now thought to be a non-opioid receptor. Mutations in this gene has been associated with juvenile amyotrophic lateral sclerosis 16. Alternative splicing of this gene results in transcript variants encoding distinct isoforms. [provided by RefSeq, Aug 2013]

**Protein Families:** Druggable Genome, GPCR, Transmembrane

### Product images:



Coomassie blue staining of purified SIGMAR1 protein (Cat# [TP301206]). The protein was produced from HEK293T cells transfected with SIGMAR1 cDNA clone (Cat# [RC201206]) using MegaTran 2.0 (Cat# [TT210002]).