

Product datasheet for **TP301207**

CCDC106 (NM_013301) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human coiled-coil domain containing 106 (CCDC106), 20 µg

Species: Human

Expression Host: HEK293T

**Expression cDNA Clone
or AA Sequence:** >RC201207 protein sequence
Red=Cloning site **Green**=Tags(s)

MNDRSSRRRTMKDDETFEISIPFDEAPHLDPQIFYLSLSPSRRNFEEPPEAASSALALMNSVKTQLHMALE
RNSWLQKRIEDLEEERDFLRCQLDKFISSARMEADHCRMKPGPRRMEGDSRGGAGGEASDPESAASSLS
GASEEGSASERRRQKQKGGASRRRFGKPKARERQRVKDADGVLCRYKKILGTFTQKLKSMSRAFEHHRVDR
NTVALTTPIAELLIVAPEKLAEVGFEFDPSEKERLLEYSRRCFLALDDETLKKVQALKKSKLLLPIYRFKR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 31.9 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_037433](#)

Locus ID: 29903


UniProt ID: [Q9BWC9](#)



[View online >](#)

RefSeq Size:	1993
Cytogenetics:	19q13.42
RefSeq ORF:	840
Synonyms:	HSU79303; ZNF581
Summary:	Promotes the degradation of p53/TP53 protein and inhibits its transactivity.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified CCDC106 protein (Cat# TP301207). The protein was produced from HEK293T cells transfected with CCDC106 cDNA clone (Cat# [RC201207]) using MegaTran 2.0 (Cat# [TT210002]).