

## Product datasheet for **TP301236L**

### GBP1 (NM\_002053) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human guanylate binding protein 1, interferon-inducible, 67kDa (GBP1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC201236 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MASEIHMTGPMCLIENTNGRLMANPEALKILSAITQPMVVAIVGLYRTGKSYLMNKLAKKKGFSLGST  
VQSHTKGIWMWCVPHPKPGHILVLLDTEGLGDVEKGDNDQNDWIFALAVLLSSTFVYNSIGTINQQAMD  
QLYYVTELTTHRIRSKSSPDENENEVEDSADVFVFPDFVWTLRDFSLDLEADGQPLTPDEYLTYSKLLKK  
GTSQKDETFNLPRLCIRKFFPKKCFVDRPVRHRRKLAQLEKLQDEELDPEFVQVADFCSYIFSNSKTK  
TLSGGIQVNGPRLESLVLTYNVAISSGDLPCMENAVLALAQIENSAAVQKAIHAYEQMGQKVLPTESL  
QELDLHRDSEREAIEVFIRSSFKDVDHLFQKELAAQLEKKRDDFCQKQNEASSDRCSGLLQVIFSPLEE  
EVKAGIYSKPGGYRLFVQKLQDLKKKYEEPRKGIQAEIILQTYLKSSESMTDAILQTDQTLTEKEKEIE  
VERVKAESAQASAKMLQEMQRKNEQMMEQKERSYQEHKQLTEKMENDRVQLLKEQERTLALKLQEQEQL  
LKEGFQKESRIMKNEIQDLQTKMRRRKACTIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	67.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online >](#)

<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_002044</a>
<b>Locus ID:</b>	2633
<b>UniProt ID:</b>	<a href="#">P32455</a>
<b>RefSeq Size:</b>	3050
<b>Cytogenetics:</b>	1p22.2
<b>RefSeq ORF:</b>	1776
<b>Summary:</b>	Guanylate binding protein expression is induced by interferon. Guanylate binding proteins are characterized by their ability to specifically bind guanine nucleotides (GMP, GDP, and GTP) and are distinguished from the GTP-binding proteins by the presence of 2 binding motifs rather than 3. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified GBP1 protein (Cat# [TP301236]). The protein was produced from HEK293T cells transfected with GBP1 cDNA clone (Cat# [RC201236]) using MegaTran 2.0 (Cat# [TT210002]).