

## Product datasheet for **TP301344**

### OSBPL2 (NM\_144498) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human oxysterol binding protein-like 2 (OSBPL2), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201344 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MNGEEFFDAVTGFSDNSSGEFSEANQKVTGMIDLDTSKNNRIGKTGERPSQENGIQKHRTSLPAPMFS  
RSDFSVWTILKKCVGLELSKITMPIAFNEPLSFLQRITEYMEHVYLIHRASCQPQPLERMQSVAFAVSA  
VASQWERTGKPFNPLLGETYELIREDLGFRFISEQVSHHPPISAFHSEGLNHDFLFHGSYIPKLFKFWGKS  
VEAEPRTITLELLKHNEAYTWTNPTCCVHNVIIGKLWIEQYGTVEILNHRTGHKCVLHFKPCGLFGKEL  
HKVEGHIQDKNKKKLFMIYGKWTECLWGIDPVSYESFKKQERRGDHLRKAKLDEDSGKADSDVADDVPA  
QETVQVIPGSKLLWRINTRPPNSAQMYNFTSFTVSLNELETGMEKTLPTDCRLRPDIRGMENGNMDLAS  
QEKERLEEKQREARRERAKEEAEWQTRWFYPGNPNPYTGTDPDWLYAGDYFERNFSDCPDIY

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	55 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.




[View online >](#)

RefSeq:	<a href="#">NP_653081</a>
Locus ID:	9885
UniProt ID:	<a href="#">Q9H1P3</a>
RefSeq Size:	4027
Cytogenetics:	20q13.33
RefSeq ORF:	1440
Synonyms:	DFNA67; DNFA67; ORP-2; ORP2

**Summary:** This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Most members contain an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain, although the encoded protein contains only the sterol-binding domain. In vitro studies have shown that the encoded protein can bind strongly to phosphatic acid and weakly to phosphatidylinositol 3-phosphate, but cannot bind to 25-hydroxycholesterol. The protein associates with the Golgi apparatus. Transcript variants encoding different isoforms have been described. [provided by RefSeq, Sep 2014]

### Product images:



Coomassie blue staining of purified OSBPL2 protein (Cat# TP301344). The protein was produced from HEK293T cells transfected with OSBPL2 cDNA clone (Cat# [RC201344]) using MegaTran 2.0 (Cat# [TT210002]).