

Product datasheet for **TP301356L**

SIGIRR (NM_021805) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens single immunoglobulin and toll-interleukin 1 receptor (TIR) domain (SIGIRR), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201356 protein sequence Red =Cloning site Green =Tags(s)

MPGVCDRAPDFLSPSEDQVLRPALGSSVALNCTAWVWSGPHCSLPSVQWLKDGLPLGIGGHYSLHEYSWW
KANLSEVLVSSVLGVNVTSTEYGAFTCSIQNISFSSFTLQ RAGPTSHVA AVLASLLVLLALLLAALLYV
KCRLNVLLWYQDAYGEVEINDGKLYDAYVSYSDCPEDRKFN FILKPKLERRRGYKFLDDRDL PRAEP
SADLLVNLSRCRRLIVLSDAFLSRAWCSHSFREGLCRLELTRRPIFITFEGQRRDPAHPALRLLRQHR
HLVTL LLWRPGSVTPSSDFWKEVQLALPRKVRYRPVEGDPQTQLQDDKDPMLILRGRVPEGRALDSEVDP
DPEGDLGVRGPVFGEPSAPPHTSGVSLGESRSSEVDVSDLGSRNYSARTDFYCLVSKDDM

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	45.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_068577</u>




[View online »](#)

Locus ID: 59307
UniProt ID: [Q6IA17](#)
RefSeq Size: 1657
Cytogenetics: 11p15.5
RefSeq ORF: 1230
Synonyms: IL-1R8; TIR8

Summary: Acts as a negative regulator of the Toll-like and IL-1R receptor signaling pathways. Attenuates the recruitment of receptor-proximal signaling components to the TLR4 receptor, probably through an TIR-TIR domain interaction with TLR4. Through its extracellular domain interferes with the heterodimerization of IL1R1 and IL1RAP.[UniProtKB/Swiss-Prot Function]

Protein Families: Transcription Factors, Transmembrane

Product images:



Coomassie blue staining of purified SIGIRR protein (Cat# [TP301356]). The protein was produced from HEK293T cells transfected with SIGIRR cDNA clone (Cat# [RC201356]) using MegaTran 2.0 (Cat# [TT210002]).