

Product datasheet for TP301362M

HMBS (NM_000190) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hydroxymethylbilane synthase (HMBS), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201362 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSGNGNAAATAEENSPKMRVIRVGTGRKSQLARIQTDSVVATLKASYPGLQFEIIMSTTGDKILDALS IGEKSLFTKELEHALEKNEVDLVVHSLKDLPTVLPFGFTIGAICKRENPHDAVWFHPKFVGGKTLLETLP SVVGTSSLRRAAQLQRKFPHEFRSIRGNLNLRLRLKDEQQEFSAILATAGLQRMGWHNVRVQILHPEE CMYAVGQGALGVEVRAKDQDILDVLVGLHDPETLLRCIAERAFLRHLEGCSVPVAVHTAMKDGQLYLTG GVWSLDGSDSIQETMQATIHVPAQHEDGPEDDPQLVGITARNIPRGPQLAAQNLGISLANLLLSKGA KNI LDVARQLNDAH</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	39.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000181



[View online »](#)


Locus ID: 3145
UniProt ID: [P08397](#)
RefSeq Size: 1526
Cytogenetics: 11q23.3
RefSeq ORF: 1083
Synonyms: PBG-D; PBGD; PORC; UPS

Summary: This gene encodes a member of the hydroxymethylbilane synthase superfamily. The encoded protein is the third enzyme of the heme biosynthetic pathway and catalyzes the head to tail condensation of four porphobilinogen molecules into the linear hydroxymethylbilane. Mutations in this gene are associated with the autosomal dominant disease acute intermittent porphyria. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways, Porphyrin and chlorophyll metabolism

Product images:



Coomassie blue staining of purified HMBS protein (Cat# [TP301362]). The protein was produced from HEK293T cells transfected with HMBS cDNA clone (Cat# [RC201362]) using MegaTran 2.0 (Cat# [TT210002]).