

## Product datasheet for **TP301368**

### RPL8 (NM\_033301) Human Recombinant Protein

#### Product data:

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human ribosomal protein L8 (RPL8), transcript variant 2, 20 µg

**Species:** Human

**Expression Host:** HEK293T

**Expression cDNA Clone or AA Sequence:** >RC201368 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MGRVIRGQRKGAGSVFRAHVKHRKGAARLRAVDFAERHGYIKGIVKDIHDPGRGAPLAKVFRDPYRFK  
KRTELFIAAEGIHGTGQFVYCGKKAQLNVGNVLPVGTMPGTVCCLEEKPGDRGKLARASGNYATVISHN  
PETKKTRVKLPSGSKKVISSANRAVVGWVAGGGRIDKPILKAGRAYHKYKAKRNCWPRVRGVMNPVEHP  
FGGGNHQHIGKPSTIRRDAPAGRKVGLIAARRTGRLRGTKTVQEKEN

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

**Tag:** C-Myc/DDK

**Predicted MW:** 27.8 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

**Storage:** Store at -80°C.

**Stability:** Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** [NP\\_150644](#)

**Locus ID:** 6132

**UniProt ID:** [P62917](#)



[View online »](#)

RefSeq Size: 967

Cytogenetics: 8q24.3


RefSeq ORF: 771

Synonyms: L8

**Summary:** Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L2P family of ribosomal proteins. It is located in the cytoplasm. In rat, the protein associates with the 5.8S rRNA, very likely participates in the binding of aminoacyl-tRNA, and is a constituent of the elongation factor 2-binding site at the ribosomal subunit interface. Alternatively spliced transcript variants encoding the same protein exist. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]

**Protein Pathways:** Ribosome

### Product images:



Coomassie blue staining of purified RPL8 protein (Cat# TP301368). The protein was produced from HEK293T cells transfected with RPL8 cDNA clone (Cat# [RC201368]) using MegaTran 2.0 (Cat# [TT210002]).