

Product datasheet for **TP301371L**

UBAP2L (NM_014847) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ubiquitin associated protein 2-like (UBAP2L), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201371 protein sequence Red =Cloning site Green =Tags(s)

MMTSVGTNRARGNWEQPQNQNQTQHKQRPQATAEQIRLAQMISDHNDADFEKVKQLIDITGKNQDECVI
ALHDCNGDVNRRAINVLLEGNPDTHSWEMVVGKKKGVSGQKDGQTESNEEGKENRDRDRDYSRRRGPPRR
GRGASRGREFRGQENGLDGTKSGGSPGRGTERGRRGRGRGGSGRRGGRFSAQGMGTFPADYAEPANT
DDNYGNSSGNTWNNTGHFEPDDGTSAWRTATEEWGTEDWNEDLSETKIFTASNVSSVPLPAENVITITAGQ
RIDLAVLLGKTPSTMENDSSNLDPSQAPSLAQLVFSNSKQTAISQPASGNTFSHHSMVSMGLGKGFVDVG
EAKGGSTTGSQFLEQFKTAQALAAQHSQSGSTTTSSWDMGTTQSPSLVQYDLKNPSDSAVHSPFTK
RQAFTPSSTMMEVFLQEKSPAVATSTAAPPPSSPLPSKSTSAPQMSPGSSDNQSSSPQPAHQKLKQKK
KASLTSKIPALAVEMPGSADISGLNLQFGALQFGSEPVLSDYESTPTTSASSSQAPSSLYTSTASESST
ISSNQSQESGYQSGPIQSTTYTSQNNAQGPLYEQRSTQTRYPSSISSSPQKDLTQAKNGFSSVQATQLQ
TTQSVGATGSAVKSDSPSTSSIPPLNETVSAASLLTTTNQHSSSLGGLSHSEEIPNTTTTQHSSTLSTQ
QNTLSSSTSSGRTSTSTLLHTSVESEANLHSSSSTFSTTSSTVSAPPPVSVSSSLNSGSSLGLSLGNS
TVTASTRSSVATTSGKAPPNLPVPLLPNPYIMAPGLLHAYPPQVYGYDDLQMLQTRFPLDYISIPFP
TPTTPLTGRDGLASNPYSGDLTKFGRGDASSPAPATLAQPQNQTQTHHTTQQTFLNPALPPGYSYTS
LPYYTGVPGLPSTFYGPVAVFPVAPTSSKQHGVNVSINASATPFQQPSGYGSHGYNTGVSVTSSNTGVPD
ISGSVYSKTQQSFEKQGFHSGTPAASFNLPSALGSGGPINPATAAAAYPPAPFMHILTPHQPHSQILHHH
LQQDGQTGSGQRSQTSSIPQKPQTNKSAYNSYSWGAN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	114.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_055662
Locus ID:	9898
UniProt ID:	Q14157
RefSeq Size:	4006
Cytogenetics:	1q21.3
RefSeq ORF:	3261
Synonyms:	NICE-4; NICE4
Summary:	Plays an important role in the activity of long-term repopulating hematopoietic stem cells (LT-HSCs).[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified UBAP2L protein (Cat# [TP301371]). The protein was produced from HEK293T cells transfected with UBAP2L cDNA clone (Cat# [RC201371]) using MegaTran 2.0 (Cat# [TT210002]).