

Product datasheet for TP301380L

FAM62A (ESYT1) (NM_015292) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 62 (C2 domain containing), member A (FAM62A), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201380 protein sequence Red=Cloning site Green=Tags(s)

MERSPGEGPSPSPMDQPSAPSDPTDQPPAAHAKPDPGSGGQPAGPGAAGEALAVLTSFGRRLLVLIPVYL
AGAVGLSVGFVLFGLALYLGWRRVRDEKERSLRAARQLLDDEEQLTAKTLYMSHRELPAWVSFPDVEKAE
WLNKIVAQVWPFLGQYMEKLLAETVAPAVRGSNPHLQFTFTFTRVELGEKPLRIIGVKVHPGQRKEQILLD
LNISYVGDVQIDVEVKKYFKAGVKGMQLHGVLRVILEPLIGDLFPVGVAVSMFFIRRPDLINWTGMTNL
LDIPGLSSLSDTMIMDSIAAFLVLPNRLVPLVPLDQDVAQLRSPLPRGIIRIHLLAARGLSKDKYVKG
LIEGKSDPYALVRLGTQTFCSRVIDEELNPQWGETYEVMMVEVPGQEIEVEVFDKDPDKDDFLGRMKLDV
GKVLQASVLDDWFPLQGGQGQVHLRLEWLSLLSDAEKLEQVLQWNWGVSSRPDPPSAAILVVYLDRAQDL
PLKKGNKEPNPMVQLSIQDVTQESKAVYSTNCPVWEEAFRFFLQDPQSQELDVQVKDDSRALTGALTLP
LARLLTAPELILDQWFQLSSSGPNSRLYMKLVMRILYLDSEICFPTVPGCPGAWDVDSNPQRGSSVDA
PPRPCHTTPDSQFGTEHVLRIHVLEAQDLIAKDRFLGGLVKGKSDPYVKLKLGRSFRSHVREDLNPRW
NEVFEVIVTSVPGQELEVEVFDKDLKDDFLGRCKVRLTTLVNSGFLDEWLTLEDVPSGRHLHLRLERLTP
RPTAAELEEVLQVNSLIQTQKSAELAAALLSIYMERAEPLRKGTKHLSPYATLTVGDSSSHKTKTISQT
SAPVWDESASFLIRKPHTESLELQVRGEGTGVLGSLPLSELLVADQLCLDRWFTLSSGQGQVLLRAQL
GILVSQHSGVEAHSHSYSHSSSSLSEPELSGGPPHITSSAPELRQLTHVDSPLEAPAGPLGQVKLTLW
YYSEERKLVSIHGCRLRQNGRDPDPYVSLLLLPDKNRGTRKRTSQQKRTLSPEFNERFEWELPLDEA
QRRKLDVSVKSNSSFMSRERELLGKVQLDLAETDLSQGVARWYDLMDNKKDKGSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	122.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_056107
Locus ID:	23344
UniProt ID:	Q9BSJ8 , A0A024RB17
RefSeq Size:	4255
Cytogenetics:	12q13.2
RefSeq ORF:	3312
Synonyms:	FAM62A; MBC2
Summary:	Binds glycerophospholipids in a barrel-like domain and may play a role in cellular lipid transport (By similarity). Binds calcium (via the C2 domains) and translocates to sites of contact between the endoplasmic reticulum and the cell membrane in response to increased cytosolic calcium levels. Helps tether the endoplasmic reticulum to the cell membrane and promotes the formation of appositions between the endoplasmic reticulum and the cell membrane. [UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transmembrane

Product images:



Coomassie blue staining of purified ESYT1 protein (Cat# [TP301380]). The protein was produced from HEK293T cells transfected with ESYT1 cDNA clone (Cat# [RC201380]) using MegaTran 2.0 (Cat# [TT210002]).