

Product datasheet for **TP301390L**

GNRPX (PLEKHJ1) (NM_018049) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human pleckstrin homology domain containing, family J member 1 (PLEKHJ1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201390 protein sequence Red =Cloning site Green =Tags(s)
	MRYNEKELQALSRQPAEMAAELGMRGPKKGSVLKRRRLVCLVNFVFLFYFRTDEAEPVGGALLLRCRVVREE PGTFSISFIEDPERKYHFECSSQQEWMEALRRASYEFMRRSLIFRNEIRKVTGKDPLEQFGISEEA RFQLSGLQA
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	17.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060519
Locus ID:	55111
UniProt ID:	Q9NW61



[View online »](#)

RefSeq Size:	1215
Cytogenetics:	19p13.3
RefSeq ORF:	447
Synonyms:	GNRPX

Product images:

Coomassie blue staining of purified PLEKHJ1 protein (Cat# [TP301390]). The protein was produced from HEK293T cells transfected with PLEKHJ1 cDNA clone (Cat# [RC201390]) using MegaTran 2.0 (Cat# [TT210002]).