

## Product datasheet for **TP301393L**

### ELOB (NM\_007108) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human transcription elongation factor B (SIII), polypeptide 2 (18kDa, elongin B) (TCEB2), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201393 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MDVFLMIRRHKTTIFTDAKESSTVFELKRIVEGILKRPPDEQRLYKDDQLLDDGKTLGECGFTSQTPARQ APATVGLAFRADDTFEALCIEPFSSPPELPDVMKPKQDSGSSANEQAVQ
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	13 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_009039</a>
Locus ID:	6923
UniProt ID:	<a href="#">Q15370</a>
RefSeq Size:	1009



[View online »](#)

Cytogenetics: 16p13.3

RefSeq ORF: 354


Synonyms: SIII; TCEB2

**Summary:** This gene encodes the protein elongin B, which is a subunit of the transcription factor B (SIII) complex. The SIII complex is composed of elongins A/A2, B and C. It activates elongation by RNA polymerase II by suppressing transient pausing of the polymerase at many sites within transcription units. Elongin A functions as the transcriptionally active component of the SIII complex, whereas elongins B and C are regulatory subunits. Elongin A2 is specifically expressed in the testis, and capable of forming a stable complex with elongins B and C. The von Hippel-Lindau tumor suppressor protein binds to elongins B and C, and thereby inhibits transcription elongation. Two alternatively spliced transcript variants encoding different isoforms have been described for this gene. Pseudogenes have been identified on chromosomes 11 and 13. [provided by RefSeq, Aug 2008]

**Protein Families:** Druggable Genome, Transcription Factors

**Protein Pathways:** Pathways in cancer, Renal cell carcinoma, Ubiquitin mediated proteolysis

### Product images:



Coomassie blue staining of purified ELOB protein (Cat# [TP301393]). The protein was produced from HEK293T cells transfected with ELOB cDNA clone (Cat# [RC201393]) using MegaTran 2.0 (Cat# [TT210002]).