

Product datasheet for TP301406M

XAB1 (GPN1) (NM_007266) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human GPN-loop GTPase 1 (GPN1), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201406 representing NM_007266 Red =Cloning site Green =Tags(s)

MAASAAAAELQASGGPRHPVCLLVLMAGSGKTTFVQRLTGHLHAQGTPPYVINLDPVHEVPFPANIDI
RDTVKYKEVMKQYGLGPNGGIVTSLNLFATRFDQVMKFIEKAQNMSKYVLIDTPGQIEVFTWSASGTIIT
EALASSFPTVVIYVMDTSRSTNPVTFMSNMLYACSIYKTKLPFIVVMNKTDIIDHSFAVEWMQDFAFQ
DALNQETTYVSNLTRSMLVLDEFYSSLRVVGVS AVLGTGLDEL FVQVTSAAEEYEREYRPEYERLKKSL
ANAESQQQREQLERLRKDMGVSVALDAGTAKDSLSPVLHPSDLILTRGTLDEEDEEADSDTDDIDHRVTEE
SHEEPAFQNFMQESMAQYWKRNNK


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	41.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_009197
Locus ID:	11321
UniProt ID:	Q9HCN4
RefSeq Size:	1829
Cytogenetics:	2p23.3
RefSeq ORF:	1122
Synonyms:	ATPBD1A; MBDIN; NTPBP; RPAP4; XAB1
Summary:	This gene encodes a guanosine triphosphatase enzyme. The encoded protein may play a role in DNA repair and may function in activation of transcription. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]

Product images:



Coomassie blue staining of purified GPN1 protein (Cat# [TP301406]). The protein was produced from HEK293T cells transfected with GPN1 cDNA clone (Cat# [RC201406]) using MegaTran 2.0 (Cat# [TT210002]).