

# Product datasheet for TP301464M

## RAB10 (NM\_016131) Human Recombinant Protein

### **Product data:**

#### **Product Type: Recombinant Proteins** Recombinant protein of human RAB10, member RAS oncogene family (RAB10), 100 µg **Description:** Species: Human HEK293T **Expression Host:** Expression cDNA Clone >RC201464 protein sequence or AA Sequence: Red=Cloning site Green=Tags(s) MAKKTYDLLFKLLLIGDSGVGKTCVLFRFSDDAFNTTFISTIGIDFKIKTVELQGKKIKLQIWDTAGQER FHTITTSYYRGAMGIMLVYDITNGKSFENISKWLRNIDEHANEDVERMLLGNKCDMDDKRVVPKGKGEQI AREHGIRFFETSAKANINIEKAFLTLAEDILRKTPVKEPNSENVDISSGGGVTGWKSKCC TRTRPLEQKLISEEDLAANDILDYKDDDDKV C-Myc/DDK Tag: Predicted MW: 22.4 kDa **Concentration:** >0.05 µg/µL as determined by microplate BCA method **Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol **Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. Storage: Store at -80°C. Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. **RefSeq:** NP 057215 Locus ID: 10890 **UniProt ID:** P61026 3785 **RefSeq Size:**



View online »


This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

	RAB10 (NM_016131) Human Recombinant Protein – TP301464M
Cytogenetics:	2p23.3
RefSeq ORF:	600
Summary:	RAB10 belongs to the RAS (see HRAS; MIM 190020) superfamily of small GTPases. RAB proteins localize to exocytic and endocytic compartments and regulate intracellular vesicle trafficking (Bao et al., 1998 [PubMed 9918381]).[supplied by OMIM, Mar 2009]
Protein Families	: Druggable Genome

## **Product images:**



Coomassie blue staining of purified RAB10 protein (Cat# [TP301464]). The protein was produced from HEK293T cells transfected with RAB10 cDNA clone (Cat# [RC201464]) using MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US