

Product datasheet for **TP301466L**

POLR1D (NM_015972) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human polymerase (RNA) I polypeptide D, 16kDa (POLR1D), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC201466 protein sequence Red =Cloning site Green =Tags(s)

MEEDQELERKISGLKTSMAEGERKTALEMVQAAGTDRHCVTFVLHEEDHTLGNSLRYMIMKNPEVEFCGY
TTTHPSESKINLRIQTRGTLPAVEPFQRGLNELMNVQCQHVLDKFEASIKDYKDQKASRNESTF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	15.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057056
Locus ID:	51082
UniProt ID:	P0DPB6
RefSeq Size:	817



[View online »](#)

Cytogenetics:	13q12.2
RefSeq ORF:	399
Synonyms:	AC19; POLR1C; RPA9; RPA16; RPAC2; RPC16; RPO1-3; TCS2
Summary:	The protein encoded by this gene is a component of the RNA polymerase I and RNA polymerase III complexes, which function in the synthesis of ribosomal RNA precursors and small RNAs, respectively. Mutations in this gene are a cause of Treacher Collins syndrome (TCS), a craniofacial development disorder. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2011]
Protein Families:	Stem cell - Pluripotency, Transcription Factors
Protein Pathways:	Cytosolic DNA-sensing pathway, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase

Product images:



Coomassie blue staining of purified POLR1D protein (Cat# [TP301466]). The protein was produced from HEK293T cells transfected with POLR1D cDNA clone (Cat# [RC201466]) using MegaTran 2.0 (Cat# [TT210002]).